




HYDRO-GAZ-MED
MEDICAL EQUIPMENT MANUFACTURER
Willowa 40 St., 05-205 Dobczyn, Poland
☎ +48 22 787 65 60 / +48 606 348 532 📠 +48 22 487 98 53
💻 www.gazmed.pl / ✉ gazmed@gazmed.pl

MEDICAL GAS FLOW MONITOR

according to EN ISO 7396-1:2016



Instruction manual

CE

Table of contents

| | |
|--|---|
| 1. General information and specification | 3 |
| 1.1. Specification | 3 |
| 2. Function | 3 |
| 3. Construction | 4 |
| 4. Screen content | 4 |
| 4.1. Main screen | 4 |
| 4.2. History | 4 |
| 5. Date & clock settings | 5 |
| 6. Transmitting alarm signals | 5 |
| 7. Maintenance | 6 |
| 8. Parts list..... | 6 |
| 9. Authorized service | 6 |
| 10. Communication cable scheme | 6 |
| 11. Wiring diagram | 7 |



IMPORTANT!
**PRIOR TO CONNECTING TO A POWER SOURCE ENSURE THAT FLOW
SENSOR IS CONNECTED TO DEVICE!**

1. General information and specification

Device is designed to measure current flow, pressure and total consumption of medical gas.

1.1. Specification

| Standard | HGM | |
|----------------------------|-----------------------|---|
| Input | Flow sensor | Thermal mass sensor FS7 |
| Output | MODBUS RTU (optional) | 9600 BAUD, 8 bits, 2 bits stop, no parity control |
| Measuring range | Flow | 0-300 l/min |
| | Pressure | 0-16 bar |
| Inlet pressure | Max. 16 bar | |
| Air inlet | 22 mm copper pipe | |
| Connectors | Transmission | RJ-45 or STL-1550/4-3.5 |
| | Inlet/outlet | STL-1550/4-3.5 MC-1.5/2-5,08 |
| Power supply | 230V AC | |
| Power consumption | ~25mA | Max. 200mA |
| Housing | PET | |
| Working temperature | 5 – 50°C | |
| Storing temperature | -20 – 60°C | |
| Dimensions | Length | 190 mm |
| | Width | 145 mm |
| | Height | 80 mm |
| Weight | ~2 kg | |

2. Function

Gas passes through flow sensor and its speed is measured and displayed as flow value. Flow is measured proportionally and displayed in liters per minute.

Device can memorize consumption of last 6 months. Total consumption can be displayed in cubic meters or in liters, this has to be determined upon order.

Precision of measurement is in range of +/- 5% within set measuring range.


3. Construction

Device consists of:

- 1) flow sensor;
- 2) alarm unit;
- 3) 12V power supply;
- 4) 230V cable;

4. Screen content

4.1. Main screen




4.2. History

History screen shows total consumption of gas for past six months. Every month is measured separately. Each consumption per month is presented in cubic meters. Total consumption will count till 655 359 m³, after that it will reset itself and start from 0.

| Monthly consumption | |
|---------------------|-------------------------------|
| Months | O ₂ m ³ |
| Present | 10 |
| 1 back | 449 |
| 2 back | 854 |
| 3 back | 630 |
| 4 back | 560 |
| 5 back | 677 |


5. Date & clock settings

Date and clock can be changed according to local time if necessary. To access date and clock settings device will ask for PIN.



Enter **5287** and tap OK.

If wrong digit or code is entered, touch back arrow and re-enter PIN again, if Pin is wrong, there will be single “beep” sound.



After adjusting tap ESC to save and go to the main screen.

6. Transmitting alarm signals

Control unit is provided with RJ-45 connectors on the back (see clause 9 and 10), they allow to connect any remote alarm unit compatible with Hydro-Gaz-Med communication interface or be connected to central monitoring system installed in healthcare facility using MODBUS RTU, such converter is sold separately.

| MODBUS | DI | RI | RO |
|--------------------|----|------------------|----|
| 9600 bps | 0. | 0. serial number | 0. |
| Device address 200 | 1. | 1. Flow value | 1. |
| Base address 0 | 2. | 2. | 2. |

7. Maintenance

Device is maintenance free.

8. Parts list

- 12V DC power supply - cat. No. APV-16-12
- flow sensor - cat. No. HGM-FS7

9. Authorized service

HYDRO-GAZ-MED Sp.J.
Willowa 40 St.
05-205 Dobczyn, Poland
tel. +48 22 787 65 60
fax +48 22 487 98 53
gazmed@gazmed.pl

10. Communication cable scheme



11. Wiring diagram

