

# Central alarm unit "S9"



Central alarm unit allows to remote monitoring multiple devices at the same time. It shows conditions, alarm states from every device within the "HGM network" connected via BMS, via MODBUS RTU or via LAN/WIFI module.

Records gas consumption or pressure values as graph.

Each central alarm unit is programmed individually for our clients, depending on their needs and the network configuration.

UI is very user-friendly and intuitive to use on 7" LCD HQ touch display.

ALL devices are connected with use of FTP cables and RJ-45 plugs.

Each alarm system can be equipped with SMS module, which allows to receive text message with alerts directly on mobile phone right after alarm appears.

Our SMS module can send messages up to 3 different numbers at the same time.