


Reducing sets for compressed medical gases



Function

To reduce pressure of compressed medical gases to a level required by user.

Technical data:

- ▶ 2 parallel pressure reducers;
- ▶ brass ball valves with chrome body and stainless steel ball;
- ▶ pressure sensor or 2 pressure switches;
- ▶ max inlet pressure 16 bar (232 psi);
- ▶ outlet pressure 1-10 bar (14-145 psi);
- ▶ max flow 50m³ or 200m³
- ▶ ambient temperature -20°C do +40°C;
- ▶ safety valve to 7 bar (101 psi) or 11 bar (159 psi);
- ▶ standard outlet pressure 5 (72 psi) and 8 bar (116 psi), others on demand.

Connections:

- ▶ copper pipe size 15 or 22 mm;